

Name: \_\_\_\_\_

Date: \_\_\_\_\_

# Microscope

## Across

**4.** Holds the high and low power objectives, can be rotated to change magnification.

**7.** Holds the slide in place.

**9.** Is used to regulate the amount of light on the specimen.

**10.** Supports the entire microscope.


**11.** Hollow tube that connects the eyepiece to the objective lenses.

**12.** Moves the stage slightly to sharpen the image.

## Down

**1.** Moves the stage up and down for focusing.

**2.** You look through to see the specimen.



**3.** Magnification of objects ranges from 10x to 40x.

**5.** Supports the slide being viewed.

**6.** Projects light upwards through the diaphragm, the specimen, and the lenses.

**8.** Is used to support the microscope when carried.

