

Name: _____

Date: _____

Microscope

Across

4. Holds the high and low power objectives, can be rotated to change magnification.

7. Holds the slide in place.

9. Is used to regulate the amount of light on the specimen.

10. Supports the entire microscope.


11. Hollow tube that connects the eyepiece to the objective lenses.

12. Moves the stage slightly to sharpen the image.

Down

1. Moves the stage up and down for focusing.

2. You look through to see the specimen.



3. Magnification of objects ranges from 10x to 40x.

5. Supports the slide being viewed.

6. Projects light upwards through the diaphragm, the specimen, and the lenses.

8. Is used to support the microscope when carried.

